

# What's New in HealthKit

Changes in iOS 9

Session 203

Shannon Tan iOS Software Engineer

Allan Shortlidge iOS Software Engineer

# Roadmap

## What's new in HealthKit

Overview

Unit preferences

New data types

Source revisions and devices

Deleted sample queries


Workout sessions

WatchKit demo









# Data Types and Units

HealthKit overview

Steps

Nike Fuel

Heart Rate

Vitamin C

Calories

Cycling Distance

Body Mass

Body Temperature

Oxygen Saturation

BMI

Blood Pressure

Vitamin B6

Potassium

Flights Climbed

Vitamin B12

BAC

Body Fat Percentage

Height

Inhaler Usage

Vitamin A

Respiratory Rate

Vitamin D

Blood Glucose

# Adding Data Types

HealthKit overview



# Adding Data Types

HealthKit overview

Existing hardware

# Adding Data Types

HealthKit overview

Existing hardware

Existing software

# Adding Data Types


HealthKit overview

Existing hardware

Existing software

Customer and developer feedback

# HKObjectType







# For More Information



# Unit Preferences

New in iOS 8.2



# Unit Preferences

```
extension HKHealthStore {  
    func preferredUnitsForQuantityTypes(quantityTypes: Set<HKQuantityType>,  
        completion: ([HKQuantityType : HKUnit], NSError?) -> Void)  
}
```

```
let HKUserPreferencesDidChangeNotification: String
```



# New Data Types

# Water Intake

## New data types

Many apps track water intake

HKQuantityTypeIdentifierDietaryWater

Volume, cumulative



# UV Exposure

New data types

HKQuantityTypeIdentifierUVExposure

Unit: Scalar (UV index)

# Fitzpatrick Skin Type

New data types

HKCharacteristicTypeIdentifierFitzpatrickSkinType

Types I-VI, or not set



I/II



III



IV



V



VI



# Fitzpatrick Skin Type

New data types

```
class HKHealthStore : NSObject {  
    ...  
    func fitzpatrickSkinTypeWithError() throws -> HKFitzpatrickSkinTypeObject  
}
```

```
class HKFitzpatrickSkinTypeObject : NSObject {  
    var skinType: HKFitzpatrickSkinType { get }  
}
```

# Reproductive Health

New data types

# Reproductive Health

## New data types

Basal body temperature

Cervical mucus quality

Ovulation

Menstrual flow

Vaginal spotting

Sexual activity





# Basal Body Temperature

New data types

```
// HKTypeIdentifiers
```

```
let HKQuantityTypeIdentifierBasalBodyTemperature: String
```





# Cervical Mucus Quality

New data types

```
// HKTypeIdentifiers
let HKCategoryTypeIdentifierCervicalMucusQuality: String

// HKDefines
enum HKCategoryValueCervicalMucusQuality : Int {
    case Dry
    case Sticky
    case Creamy
    case Watery
    case EggWhite
}
```



# Ovulation Test Result

New data types

```
// HKTypeIdentifiers
let HKCategoryTypeIdentifierOvulationTestResult: String

// HKDefines
enum HKCategoryValueOvulationTestResult : Int {
    case Negative
    case Positive
}
```



# Menstruation

## New data types

```
// HKTypeIdentifiers
let HKCategoryTypeIdentifierMenstrualFlow: String

// HKDefines
enum HKCategoryValueMenstrualFlow : Int {
    case Unspecified
    case Light
    case Medium
    case Heavy
}
```



# Menstruation: Metadata

New data types

```
// HKMetadata
// The expected value is a Boolean
// Required
let HKMetadataKeyMenstrualCycleStart: String
```

# Metadata Example

## New data types

```
let healthStore: HKHealthStore = HKHealthStore()

let dict = [HKMetadataKeyMenstrualCycleStart: true]

let type = HKCategoryType.categoryTypeForIdentifier(
    HKCategoryTypeIdentifierMenstrualFlow)!
let value = HKCategoryValueMenstrualFlow.Unspecified.rawValue
let date = NSDate()

let sample = HKCategorySample(type: type, value: value, startDate: date,
    endDate: date, metadata: dict)

healthStore.saveObject(sample) { ... }
```

# Metadata Example

## New data types

```
let healthStore: HKHealthStore = HKHealthStore()
```

```
let dict = [HKMetadataKeyMenstrualCycleStart: true]
```

```
let type = HKCategoryType.categoryTypeForIdentifier(  
    HKCategoryTypeIdentifierMenstrualFlow)
```

```
let value = HKCategoryValueMenstrualFlow.Unspecified.rawValue
```

```
let date = NSDate()
```

```
let sample = HKCategorySample(type: type, value: value, startDate: date,  
    endDate: date, metadata: dict)
```

```
healthStore.saveObject(sample) { ... }
```

# Metadata Example

## New data types

```
let healthStore: HKHealthStore = HKHealthStore()

let dict = [HKMetadataKeyMenstrualCycleStart: true]

let type = HKCategoryType.categoryTypeForIdentifier(
    HKCategoryTypeIdentifierMenstrualFlow)!
let value = HKCategoryValueMenstrualFlow.Unspecified.rawValue
let date = NSDate()

let sample = HKCategorySample(type: type, value: value, startDate: date,
    endDate: date, metadata: dict)

healthStore.saveObject(sample) { ... }
```



# Single-Value Category Samples

## New data types

Next two types will be category samples

New category sample value: `HKCategoryValueNotApplicable`

Used for category samples that have dates, metadata, and other attributes, but do not have multiple values

# Vaginal Spotting

New data types

```
// HKTypeIdentifiers
```

```
let HKCategoryTypeIdentifierVaginalSpotting: String
```



# Sexual Activity

New data types

```
// HKTypeIdentifiers
let HKCategoryTypeIdentifierSexualActivity: String

// HKMetadata
// The expected value is a Boolean
let HKMetadataKeySexualActivityProtectionUsed: String
```



# Summary

New data types

Water intake

UV exposure/Fitzpatrick Skin Type

Six reproductive health types

# Data Sources

















# Sources and Device Metadata

## Data sources

```
// HKObject
class HKObject {
    var source: HKSource { get }
    var metadata: [String : AnyObject]? { get }

}

// HKMetadata
let HKMetadataKeyDeviceName: String
let HKMetadataKeyDeviceManufacturerName: String
...
```



# Sources and Device Metadata

## Data sources

```
// HKObject
class HKObject {
    var source: HKSource { get }
    var metadata: [String : AnyObject]? { get }
}

// HKMetadata
let HKMetadataKeyDeviceName: String
let HKMetadataKeyDeviceManufacturerName: String
...
```

# Sources and Device Metadata

## Data sources

```
// HKObject
class HKObject {
    var source: HKSource { get }
    var metadata: [String : AnyObject]? { get }
}

// HKMetadata
let HKMetadataKeyDeviceName: String
let HKMetadataKeyDeviceManufacturerName: String
...
```

# Object Source Revisions and Devices

Data sources

// HKObject

```
class HKObject {
```

```
    var source: HKSource { get }
```

```
    var metadata: [String : AnyObject]? { get }
```

```
    var sourceRevision: HKSourceRevision { get }
```

```
    var device: HKDevice? { get }
```

```
}
```

// HKMetadata

```
let HKMetadataKeyDeviceName: String
```

```
let HKMetadataKeyDeviceManufacturerName: String
```

```
...
```

# Source Revisions

Data sources

```
class HKSourceRevision {  
    var source: HKSource { get }  
    var version: String? { get } // kCFBundleVersionKey  
}
```

# Devices

## Data sources

```
class HKDevice {  
    var name: String? { get }  
    var manufacturer: String? { get }  
    var model: String? { get }  
    var hardwareVersion: String? { get }  
    var firmwareVersion: String? { get }  
    var softwareVersion: String? { get }  
    var localIdentifier: String? { get }  
    var UUIDeviceIdentifier: String? { get }  
  
}
```

# Devices

## Data sources

```
class HKDevice {  
    var name: String? { get }  
    var manufacturer: String? { get }  
    var model: String? { get }  
    var hardwareVersion: String? { get }  
    var firmwareVersion: String? { get }  
    var softwareVersion: String? { get }  
    var localIdentifier: String? { get }  
    var UUIDeviceIdentifier: String? { get }  
  
    class func localDevice() -> HKDevice  
}
```



# Saving Objects with Devices

## Data sources

```
let device = HKDevice(name: "Scale", manufacturer: ... )
```

```
let sample = HKQuantitySample(type: quantityType,  
    quantity: quantity,  
    startDate: startDate,  
    endDate: endDate,  
    device: device,  
    metadata: nil)
```

```
healthStore.saveObject(sample) { ... }
```

# Saving Objects with Devices

Data sources

```
let device = HKDevice(name: "Scale", manufacturer: ... )
```

```
let sample = HKQuantitySample(type: quantityType,  
    quantity: quantity,  
    startDate: startDate,  
    endDate: endDate,  
    device: device,  
    metadata: nil)
```

```
healthStore.saveObject(sample) { ... }
```

# Saving Objects with Devices

Data sources

```
let device = HKDevice(name: "Scale", manufacturer: ... )
```

```
let sample = HKQuantitySample(type: quantityType,  
    quantity: quantity,  
    startDate: startDate,  
    endDate: endDate,  
    device: device,  
    metadata: nil)
```

```
healthStore.saveObject(sample) { ... }
```

# Saving Objects with Devices

## Data sources

```
let device = HKDevice(name: "Scale", manufacturer: ... )
```

```
let sample = HKQuantitySample(type: quantityType,  
    quantity: quantity,  
    startDate: startDate,  
    endDate: endDate,  
    device: device,  
    metadata: nil)
```

```
healthStore.saveObject(sample) { ... }
```

# Querying for Devices and Source Revisions

Data sources

```
extension HKQuery {  
    class func predicateForObjectsFromSourceRevisions(  
        sourceRevisions: Set<HKSourceRevision>) -> NSPredicate  
  
    class func predicateForObjectsFromDevices(  
        devices: Set<HKDevice>) -> NSPredicate  
}
```

# Querying for Devices and Source Revisions

## Data sources

```
let source = HKSource.defaultSource()
let revision = HKSourceRevision(source: source, version: "1.0")

let pred =
HKQuery.predicateForObjectsFromSourceRevisions([revision])

let query = HKSampleQuery(sampleType: sampleType,
    predicate: pred,
    limit: 0,
    sortDescriptors: nil) { ... }

healthStore.executeQuery(query)
```

# Querying for Devices and Source Revisions

Data sources

```
let source = HKSource.defaultSource()  
let revision = HKSourceRevision(source: source, version: "1.0")
```

```
let pred =  
HKQuery.predicateForObjectsFromSourceRevisions([revision])
```

```
let query = HKSampleQuery(sampleType: sampleType,  
    predicate: pred,  
    limit: 0,  
    sortDescriptors: nil) { ... }
```

```
healthStore.executeQuery(query)
```



# Querying for Devices and Source Revisions

## Data sources

```
let source = HKSource.defaultSource()
let revision = HKSourceRevision(source: source, version: "1.0")
```

```
let pred =
HKQuery.predicateForObjectsFromSourceRevisions([revision])
```

```
let query = HKSampleQuery(sampleType: sampleType,
    predicate: pred,
    limit: 0,
    sortDescriptors: nil) { ... }
```

```
healthStore.executeQuery(query)
```

# Querying for Devices and Source Revisions

## Data sources

```
let source = HKSource.defaultSource()  
let revision = HKSourceRevision(source: source, version: "1.0")
```

```
let pred =  
HKQuery.predicateForObjectsFromSourceRevisions([revision])
```

```
let query = HKSampleQuery(sampleType: sampleType,  
    predicate: pred,  
    limit: 0,  
    sortDescriptors: nil) { ... }  
  
healthStore.executeQuery(query)
```

Deletion

# Deleting Objects

# Deletion

```
class HKHealthStore {
    func deleteObject(object: HKObject,
        completion: (Bool, NSError?) -> Void)

}
```

# Deleting Multiple Objects

## Deletion

```
class HKHealthStore {  
    func deleteObject(object: HKObject,  
        completion: (Bool, NSError?) -> Void)  
  
    func deleteObjects(objects: [HKObject],  
        completion: (Bool, NSError?) -> Void)  
  
}
```

# Deleting Multiple Objects

## Deletion

```
class HKHealthStore {  
    func deleteObject(object: HKObject,  
        completion: (Bool, NSError?) -> Void)  
  
    func deleteObjects(objects: [HKObject],  
        completion: (Bool, NSError?) -> Void)  
  
    func deleteObjectsOfType(objectType: HKObjectType,  
        predicate: NSPredicate,  
        completion: (Bool, UInt, NSError?) -> Void)  
}
```

# Querying for Deleted Objects

## Deletion

No easy way to determine which samples were deleted with iOS 8

Difficult to synchronize another database with HealthKit



# Anchored Queries for Deleted Objects

## Deletion

```
class HKAnchoredObjectQuery {  
    init(type: HKSampleType,  
        predicate: NSPredicate?,  
        anchor: Int,  
        limit: UInt,  
        resultsHandler: (HKAnchoredObjectQuery, [HKSample]?,  
                        [HKDeletedObject]?, Int, NSError?) -> Void)  
  
}
```

# Deleted Objects

## Deletion

```
class HKDeletedObject {  
    var UUID: NSUUID { get }  
}
```

# Background Delivery of Deleted Objects

## Deletion

```
class HKHealthStore {  
    func enableBackgroundDeliveryForType(type: HKObjectType,  
        frequency: HKUpdateFrequency,  
        completion: (Bool, NSError?) -> Void)  
}
```

# Streaming Updates

## Deletion

```
class HKAnchoredObjectQuery {
    init(type: HKSampleType,
        predicate: NSPredicate?,
        anchor: Int,
        limit: UInt,
        resultsHandler: (HKAnchoredObjectQuery, [HKSample]?,
            [HKDeletedObject]?, Int, NSError?) -> Void)

    var updateHandler: ((HKAnchoredObjectQuery, [HKSample]?,
        [HKDeletedObject]?, Int, NSError?) -> Void)?
}
```

# Streaming Updates

## Deletion

```
class HKAnchoredObjectQuery {  
    init(type: HKSampleType,  
        predicate: NSPredicate?,  
        anchor: Int,  
        limit: UInt,  
        resultsHandler: (HKAnchoredObjectQuery, [HKSample]?,  
                        [HKDeletedObject]?, Int, NSError?) -> Void)  
  
    var updateHandler: ((HKAnchoredObjectQuery, [HKSample]?,  
                        [HKDeletedObject]?, Int, NSError?) -> Void)?  
}
```

HealthKit for watchOS

# Using HealthKit in Apple Watch Apps

## HealthKit on watchOS

Same APIs as iOS (limited historical data)



Access to activity and workout data

New workout APIs

Data syncs to companion device








# Privacy

## HealthKit on watchOS

WatchKit apps must request authorization

User prompted on companion device

Authorizations shared between devices




# Authorization

## HealthKit on watchOS

WatchKit apps must request authorization

User prompted on companion device

Authorizations shared between devices



# Recording Workouts

## HealthKit on watchOS

Use HKWorkoutSession while recording a workout

Specify activity type to improve accuracy

Apps stay foregrounded while in session

Only one workout session may run at a time



# Workout Sessions

HealthKit on watchOS

```
class HKWorkoutSession : NSObject {  
    var activityType: HKWorkoutActivityType { get }  
    var locationType: HKWorkoutSessionLocationType { get }  
    weak var delegate: HKWorkoutSessionDelegate?  
  
    init(activityType: HKWorkoutActivityType,  
         locationType: HKWorkoutSessionLocationType)  
}
```

# Workout Sessions

HealthKit on watchOS

```
class HKWorkoutSession : NSObject {  
    var activityType: HKWorkoutActivityType { get }  
    var locationType: HKWorkoutSessionLocationType { get }  
    weak var delegate: HKWorkoutSessionDelegate?  
  
    init(activityType: HKWorkoutActivityType,  
         locationType: HKWorkoutSessionLocationType)  
}
```

# Workout Sessions

HealthKit on watchOS

```
class HKWorkoutSession : NSObject {  
    var activityType: HKWorkoutActivityType { get }  
    var locationType: HKWorkoutSessionLocationType { get }  
    weak var delegate: HKWorkoutSessionDelegate?  
  
    init(activityType: HKWorkoutActivityType,  
         locationType: HKWorkoutSessionLocationType)  
}
```



# Workout Session Delegate

HealthKit on watchOS

```
protocol HKWorkoutSessionDelegate {  
    func workoutSession(workoutSession: HKWorkoutSession,  
        didChangeToState toState: HKWorkoutSessionState,  
        fromState: HKWorkoutSessionState,  
        date: NSDate)  
  
    func workoutSession(workoutSession: HKWorkoutSession,  
        didFailWithError error: NSError)  
}
```

# Workout Session Delegate

HealthKit on watchOS

```
protocol HKWorkoutSessionDelegate {  
    func workoutSession(workoutSession: HKWorkoutSession,  
        didChangeToState toState: HKWorkoutSessionState,  
        fromState: HKWorkoutSessionState,  
        date: NSDate)  
  
    func workoutSession(workoutSession: HKWorkoutSession,  
        didFailWithError error: NSError)  
}
```

# Workout Session Delegate

HealthKit on watchOS

```
protocol HKWorkoutSessionDelegate {  
    func workoutSession(workoutSession: HKWorkoutSession,  
        didChangeToState toState: HKWorkoutSessionState,  
        fromState: HKWorkoutSessionState,  
        date: NSDate)  
  
    func workoutSession(workoutSession: HKWorkoutSession,  
        didFailWithError error: NSError)  
}
```

# Starting and Stopping Workout Sessions

HealthKit on watchOS

```
healthStore.startWorkoutSession(workoutSession) { ... }
```

```
...
```

```
healthStore.stopWorkoutSession(workoutSession) { ... }
```

*Demo*

# Demo Recap

## HealthKit on watchOS

Requested authorization for workout data types

Started an HKWorkoutSession

Streamed samples with HKAnchoredObjectQuery

Saved HKWorkout and associated samples

# Summary

Unit preferences

New data types

Source revisions and devices

Deleted sample queries

Workout sessions

WatchKit demo



# More Information

## Documentation

<http://developer.apple.com/healthkit>

## Technical Support

Apple Developer Forums

<http://developer.apple.com/forums>

Developer Technical Support

<http://developer.apple.com/support/technical>

## General Inquiries

[healthkit@apple.com](mailto:healthkit@apple.com)

# Related Sessions and Labs

HealthKit and ResearchKit Lab	Frameworks Lab B	Wednesday 11:00AM
Building Apps with ResearchKit	Mission	Wednesday 4:30PM
Health, Fitness, and Research Get Together	Buena Vista Park	Wednesday 6:00PM
HealthKit and ResearchKit Lab	Frameworks Lab C	Thursday 11:00AM

